



# FLEXCELL® BIOFLEX® CELL SEEDERS

Optimize cell placement in each BioFlex® well!

- Confines cells during plating to the central area of a BioFlex® well (Fig. 24).
- Plating cells with a Cell Seeder allows for cells to be subjected to defined strains when using the equibiaxial Loading Stations with the Flexcell® Tension System.



**Figure 24.** The results of using a BioFlex® Cell Seeder when plating cells to confine them to the area directly above the 25 mm cylindrical loading posts (purple = cells stained with Crystal Violet stain).

#### ORDERING INFORMATION (Please contact distributor for pricing)

Catalog Number	Product/Item	Price
BFCS-1000	BioFlex Cell Seeder	
BFCS-4000	BioFlex Cell Seeder (Set of 4)	