FLEXCELL® TRANSWELL HOLDER


Holds Transwells to allow for cell migration and co-culture assays to be performed in Flexcell® culture plates.

- Analyze cell migration in response to strain.
- Run experiments with the unstrained control and the strained cells in the same culture plate well.
- Available for both 6- and 24-well Flexcell[®] culture plates.
- Compatible with standard Transwell sizes for 6-well and 24 -well plates.

ORDERING INFORMATION (Please contact distributor for pricing)		
Catalog Number	Product/Item	Price
TW6-3000	Transwell Holder for 6-well plates (Set of 4)	
TW24-3000	Transwell Holder for 24-well plates (Set of 4)	

*Transwell is a registered trademark of Corning Inc.








Figure 28. Schematic of a Transwell Holder with a Transwell supported atop the rubber membrane in the well of BioFlex® culture plate.