


CULTURE SLIPS®

**Microscope slides and cover slips surface treated for cell attachment.
Choose from any of Flexcell's culture surface treatments.**

- Available in two sizes: 75 mm × 25 mm × 1.0 mm (for Streamer or FlexFlow), 75 mm × 24 mm × 0.2 mm (for FlexFlow only)
- 25 mm × 75 mm × 1.0 mm Culture Slips are rimmed with a 1.0 mm wide Teflon border to help limit cell culture growth to the portion of the slip exposed to fluid flow.
- Low autofluorescence
- Matrix-treated to promote cell attachment.
- Covalently bonded surfaces: Amino, Collagen (Type I or IV) Elastin, ProNectin (RGD), Laminin (YIGSR)
- Delivered in sterile twin packs for one time immediate use.



ORDERING INFORMATION (call your distributor for prices)	
75 mm × 25 mm × 1.0 mm for Streamer or FlexFlow	
<i>Catalog Number</i>	<i>Product/Item</i>
CS-U	Culture Slips — Untreated
CS-A	Culture Slips — Amino
CS-C	Culture Slips — Collagen Type I
CS-C(IV)	Culture Slips — Collagen Type IV
CS-E	Culture Slips — Elastin
CS-P	Culture Slips — ProNectin
CS-L	Culture Slips — Laminin
75 mm × 24 mm × 0.2 mm for FlexFlow only	
<i>Catalog Number</i>	<i>Product/Item</i>
FFCS-U	Culture Slips — Untreated
FFCS-A	Culture Slips — Amino
FFCS-C	Culture Slips — Collagen Type I
FFCS-C(IV)	Culture Slips — Collagen Type IV
FFCS-E	Culture Slips — Elastin
FFCS-P	Culture Slips — ProNectin
FFCS-L	Culture Slips — Laminin