


# STAGEFLEXER<sup>®</sup> MEMBRANES

Flexible growth surface for the StageFlexer and FlexFlow devices

- Individually packaged 43 mm silicone rubber membranes.
- Membranes come in sterile culture dishes.
- Clear for direct viewing of cells optically.
- Same growth surfaces as the BioFlex culture plates.
- Covalently bonded surfaces: Amino, Collagen (Type I or IV) Elastin, ProNectin (RGD), Laminin (YIGSR)



ORDERING INFORMATION (call your distributor for prices)	
<i>Catalog Number</i>	<i>Product/Item</i>
<b>SFM-U</b>	StageFlexer Membrane — Untreated
<b>SFM-A</b>	StageFlexer Membrane — Amino
<b>SFM-C</b>	StageFlexer Membrane — Collagen Type I
<b>SFM-C(IV)</b>	StageFlexer Membrane — Collagen Type IV
<b>SFM-E</b>	StageFlexer Membrane — Elastin
<b>SFM-P</b>	StageFlexer Membrane — ProNectin
<b>SFM-L</b>	StageFlexer Membrane — Laminin